

Beautifully Intelligent Technology

TLC Monitoring:

- Therapeutic Logic Control guarantees consistent treatments time after time.

Temperature Monitoring:

- Real time temperature is displayed on the device, with temperature measurement response time up to 100x quicker than the competition.¹

High Patient Satisfaction:

- 99% of patients described the treatment as comfortable or relaxing.¹ With four soothing treatment tones to choose from, your patients will look forward to returning for follow-up visits.

TempSure™ *Envi*



TempSure Generator Specs & Accessories

TempSure Generator	
Size:	18"W x 22.5"D x 12"H (w/cart <46"H)
Weight:	30 lbs. (15kgs) for the Generator
Display:	10 ± 1" diagonal w/touch screen
Frequency:	4.0 MHz Monopolar
Electrical:	100-240 VAC, 50/60 Hz
Output Power:	MAX 300 watts
Face & Body Handpieces	
10mm	100-7027-110
15mm	100-7027-115
20mm	100-7027-120
25mm	100-7027-125
30mm	100-7027-130
Massage Heads	
25 mm	100-7027-055
30 mm	100-7027-056

CYNOSURE
A Hologic Company

Corporate Headquarters
5 Carlisle Road • Westford, Massachusetts 01886 USA
TELEPHONE 978 256 4200 TOLL FREE 800 886 2966
www.cynosure.com

*Through soft-tissue coagulation. ¹In-House Study. ²State laws vary with regard to delegation of RF treatments to non-physicians, and change over time. Please consult local regulations. ©2018 Hologic, Inc. All Rights Reserved. Cynosure is a registered trademark of Cynosure, Inc. TempSure is a trademark of Cynosure, Inc. Beautifully Intelligent is a trademark of Cynosure Hologic owns exclusive rights to photography. Use of photography without written permission of Hologic is prohibited. Model for illustrative purposes only. Results may vary.

CYNOSURE
A Hologic Company



TempSure™ *Envi*
Beautifully Intelligent™ RF Technology

921-7027-012 Rev.1 - 1/18

Powered by Beautifully Intelligent Technology



Introducing TempSure, Cynosure's new state-of-the-art 300W monopolar radiofrequency platform.

- Unique time & temperature monitoring (TLC) will only activate the treatment timer when the target tissue is at therapeutic temperature, delivering repeatable treatment results procedure after procedure
- Powerful 300W generator for more consistent bulk heating during body treatments
- Expandable platform that will support additional treatment modalities as your practice grows
- TempSure delivers safe, consistent, reliable RF skin tightening* treatments in any practitioner's hands²



Five HP sizes ranging from 10mm–30mm
Add-on feature for cellulite treatment

TempSure Envi face and body application reduces wrinkles and tightens skin* by non-invasively heating the deep layers of the dermis without damaging the epidermis. This encourages collagen contraction and new collagen formation.



“It doesn't start counting down the time until you get to temperature; and I think that's huge! ... And that's why it relates to results.”

—R. Saluja, Saluja Cosmetic & Laser Center

With TempSure Envi, you can give your patients an easy, gentle way to maintain tighter, younger-looking skin without invasive procedures or downtime, no matter what their skin type or the time of year. And you can give your business a way to maximize its ROI that only increases over time.

That's the beauty of TempSure Envi.

The Proof is in the Patients



Before, and one month after 2nd treatment. Courtesy of S. Doherty, MD.

For more information please visit www.cynosure.com/TempSure